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D

Multisil-Mask weich

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GB

Multisil-Mask soft

Before using the product, please read these instructions for use carefully.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product

Multisil-Mask soft is an addition-linking silicone which is injected from the double-mix cartridge directly into the impression. This way reproduction of the gingiva on the working model is achieved.

2. Indication range

Multisil-Mask soft is used for the reconstruction of the gingival margin after exposing the preparation margin on the sawcut model. It serves as an aid for gap-free modelling resp. layering of fixed dentures towards the gingival margin.

3. Contraindication

Based on the current state of knowledge there are no known contraindications.

4. Hazard and safety information

Multisil sealing agent:
After skin contact, wash with copious amount of water and soap. Provide fresh air to persons who have inhaled vapours. After eye contact, rinse with copious amount of water.

5. Interaction with other materials

The highly reactive catalyst that is contained in Multisil-Mask soft just like in all addition-linking silicones reacts highly sensitive to certain substances („catalyst poisons“). The most important catalyst poisons are organic sulphur compounds, amines and amides and compounds containing heavy metals such as lead, tin, cadmium or mercury. Contact of Multisil-Mask soft which has not hardened with such compounds (e.g. with latex gloves, partly also with condensation-linking silicones) results in damage to the setting behaviour (inhibition, retarded hardening, formation of a grease layer) and must be avoided.

6. Storage and durability

Store Multisil-Mask soft at room temperature (not above 28 °C). Protect against sunlight. Multisil-Mask soft has a shelf life of 2 years if it is stored in the original, sealed container.

7. Processing

7.1 Preparation of the model

The arch is prepared and a base is produced in the usual way. The unsawed gingival pattern is duplicated with addition – or condensation-linking kneading silicone and in this way the situation for the preparation of the gingival mask is preserved. Saw the arch and expose the die preparation. Saw cuts are filled with wax to avoid the penetration of Mask silicone. Remove excess wax to ensure correct repositioning of the silicone matrix. Prepare the silicone matrix for the application of Multisil-Mask soft. Drill inlet and outlet holes into the matrix. Do not drill in the direct vicinity of the gingival margin! Insulate the silicone matrix with Multisil Sep and reposition the matrix on the model. When repositioning, it must be ensured that separating liquid is applied into the inlet and outlet holes to avoid uncontrolled breaking of the injection channels from the actual gingival mask. Attach the matrix with sticky wax (REF 510 0040 O).

7.2 Injecting Multisil-Mask

Insert the cartridge into the dispensing device. Remove the cap of the cartridge and replace it with a mixing cannula. Dispensing and applying are simply carried out by pressing the lever of the dispensing device. Manual mixing of the components is not required. Since – in the case of cartridge materials – it is not possible to spread a 1:1 mixture of both components from the beginning, discard a small quantity (size of a hazelnut) of the material that is dispensed first. Hold the tip of the mixing cannula to an inlet hole and inject the mask material until some of the material flows out of the outlet holes. At room temperature Multisil-Mask soft can be processed approx. 1 minute. Once the application is completed, let Multisil-Mask soft harden under atmospheric pressure. Deflask the gingival mask 6 to 7 minutes after the beginning of the injection process.

7.3 Finishing the gingival mask

After lifting off the matrix, remove excess material from the gingival mask with surgical scissors or a silicone bur. For this purpose the gingival mask must be carefully lifted off the dies. Apply a thin coat of Multisil sealing agent to seal the surface. Please observe the instructions for use of Multisil sealing agent.

8. Technical data

Processing time incl. mixing time:	approx. 1 minute
Time until material can be deflasked:	approx. 6 to 7 minutes
Final hardness:	approx. 35 Shore A
Polymerization shrinkage:	< 0.1 %
Compression set:	< 1.0 %

9. Symbols

	Can be used until YYYY-MM
	Lot number

10. Additional information

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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